


ColloSkin-M

Sterile Wet Collagen Sheet in Meshed Form



Introduction

Colloskin-M is a meshed sterile collagen in sheet form. Colloskin-M's use has been widely accepted biological wound dressing especially for burns. It is highly effective for patients including those with HIV and Hepatitis infections.

Available Sizes

- 5 X 5 cms
- 10 X 10 cms
- 10 X 25 cms
- 15 X 30 cms

Advantages And Directions Of Use

Similar To Colloskin

Indications

- Idle For Exudate/infected burns and wounds for easy drainage of fluid.
- Superficial & Deep Dermal Burns.
- For 3rd Degree Burns as a temporary cover.
- Traumatic loss of skin cover.
- Temporary wound cover in major open fracture wound(s).
- Skin donor sites

60% Burn Injury. 1st Day
dressing with Colloskin.
Single application



3rd Day. Colloskin
dried up and stuck to
the wound area.



21st Day. The Wound has
healed completely and the
colloskin peeled off totally.

