

MINI-COG

Radio Frequency Electro surgical Cautery Unit

- Scars
- Hybrid skin resurfacing
- skin rejuvenation
- Fine wrinkled
- Enlarged pores
- Stretch mark
- Melasma


Accessories



Single Foot Switch



Reusable Foot Switch Pencil



1. Needle tip(Micro)
2. Wire Needle tip
3. Needle tip(45°)

3ea Tips



Bipolar Forcep



Bipolar Forcep Cable

A skin, plastic surgery operation



Nevus and Freckle Removal
Vaporization by using neeble electrode
hardly leaves a scar because the
surrounding tissue is not damaged



Ophthalmology



Power Consumption 150VA
Power Source 100-250V / 50-60Hz
Protection Class/Type Class I / B
Monopolar 35W (300Ω)

Bipolar 35W (300Ω)
Weight 2.2kg
Dimension (W×D×H) 195×165×66(mm)

itc21.co.kr

310-7, Techno2-ro, Yuseong-gu, Deajeon, Korea. / Tel : +82-42-626-0851 / Fax : +82-624-0129 / E-mail: itc@itc21.co.kr
대전광역시 유성구 테크노 2로 310-7 / Tel : 042-625-0852(1) / Fax : 042-624-0129 / Email : itc21@itc21.co.kr

